

# FILT'RAY<sup>2G</sup>

## In-line filters

AGAINST WATERBORNE  
CONTAMINANTS

SINGLE  
USE  
↓  
**62**  
DAYS

SINGLE  
USE  
↓  
**31**  
DAYS

**QUICK LEGIONELLA  
& PSEUDOMONAS CONTROL**



### INNOVATION

→ + attractive,  
+ ergonomic,  
+ functional  
for discerning users

→ A range of rigorously  
consistent filters



#### Unparalleled performance

Bubl'air Wash™ Patented  
Self-cleaning membrane  
UltraResistant (UR) high flow rate  
Tubular membrane

#### Attractive design

Compact in size  
Smooth, rounded edges

#### Ecological footprint


Recyclable materials  
"Economical" water flow rate  
Manufactured in France

## Technological Innovation - FILT'RAY<sup>26</sup> Filters

FILT'RAY<sup>26</sup> Point-of-Use filters provide a barrier against waterborne microorganisms (aerobic flora that grow at 22°C to 36°C, *Pseudomonas aeruginosa*, *Legionella pneumophila* etc.).

31 & 62 days FILTRAY<sup>26</sup> Point-of-Use filters are intended for rinsing medical devices such as flexible endoscopes. These filters are tools for environmental management and they are used to protect patients and/or residents against waterborne microorganisms.

Filters from the FILT'RAY<sup>26</sup> range are CE Class 1 medical devices that have received a Certificate of Sanitary Compliance and WRAS\*, KTW\*, and NSF\* validations. (\* under validation)



FILT'RAY<sup>26</sup> filters have tolerance to early blockage due to a **patented membrane selfcleaning system Publ'air Wash™** and provide high quality filtration and a comfortable flow rate for users.

The **Ultraresistant tubular membrane** provides sterilizing-grade microfiltration and preserves the water's sensory qualities.

### General information

- Instant microbial challenge at D0, 31 & 62 days
  - with ***Brevundimonas diminuta*** at 10<sup>7</sup> UFC/cm<sup>2</sup> of effective filtration area ; > log7 retention. The membrane is therefore considered sterilizing-grade.
  - with ***Legionella pneumophila* and *Pseudomonas aeruginosa*** at 10<sup>7</sup> UFC/cm<sup>2</sup> of effective filtration area ; > log7 retention
- Quality control : continuous complete monitoring of the filtration membrane
- Materials selected in accordance with European and US Pharmacopeia and FDA requirements
- Integrated bacteriostatic agent for a more than 99% reduction in bacterial counts
- Sterile filters sterilized by gamma irradiation in compliance with standard NF EN ISO 11137

## innovation

### Unparalleled performance

- Publ'air Wash™ Patented Self-cleaning membrane
- UltraResistant (UR) high flow rate tubular membrane
- Long-lasting filter
- 3 x larger microfiltration surface area

### Attractive appearance for real-world use

- Compact in size
- Soft, rounded edges
- Protective covering
- Easy to install

### Ecological footprint

- Recyclable materials
- "Economic" flow rate

### Component disposal policy

- Can be disposed of with normal household waste.



Characteristics	FILT'RAY <sup>26</sup> Shower head
Device materials	MBS
Membrane type	tubular microfiltration
Membrane materials	Polyethylene
Pore size	0.1 µm
Filtration flow rate	6.3 L/min at 3 bar – 2.2 L/min at 1 bar
Range of use	5 bars at 60°C
Resistance to thermal shock	70°C for 30 minutes
Resistance to chemical shock	100 ppm active chlorine during a cumulative total of 1 hour (31 Days filter) 100 ppm active chlorine during a cumulative total of 2 hours (62 Days filter)

Ref.	Product	Units
AT030733	31 Days - FILT'RAY <sup>26</sup> in-line water filter	Box of 10
AT030734	62 Days - FILT'RAY <sup>26</sup> in-line water filter	Box of 10

FILT'RAY <sup>26</sup> F Filters	31 days	62 days
Filtration membrane surface area	2800 cm <sup>2</sup>	
Membrane thickness	85 µm	110 µm

Scientific and technical validation guide for this product available upon request.

FILT'RAY<sup>26</sup> filters are manufactured in France by **aqua-tools**, 26, rue Charles-Édouard Jeanneret 78300 Poissy – France

Contact-us :  
Tél. : +33 1 39 75 02 20  
Fax : +33 1 39 75 08 28  
e-mail : [contact@aqua-tools.com](mailto:contact@aqua-tools.com)  
[www.aqua-tools.com](http://www.aqua-tools.com)

