

# FILT'RAY 1-month tap filter

Ref. AT030728

Non-sterile



## DESCRIPTION

**FILT'RAY 1-month tap filter**

**Ref. AT030728**

AQUATOOLS anti-Legionella and all-germ tap filter, shower spray.  
Device made from MBS.  
Membrane type: tubular microfiltration.  
Membrane made from polyethylene.  
Pore cut-off threshold: 0.1µm nominal and 0.2µm absolute sterilising grade.  
Filtration flow rate: 6 lpm at 3 bar and 4 lpm at 1 bar (filter only, excludes any flow rate restriction in the tap).  
Maximum upstream pressure at the point-of-use: 5 bar.  
Maximum lifespan: 1 month from installation.  
Withstands thermal shocks: 70°C for 30 minutes.  
Withstands chemical shocks: active chlorine levels up to 100 ppm for 1 hour (cumulative total).  
Non-sterile filter.  
CE marked Class I medical device.  
Packaged in boxes of 10, each filter is individually wrapped and sterile.

## ADVANTAGES



Anti-Legionella and all-germ



Lifespan: 1 month from installation



Bubl'air Wash: longer filter life



Rapid installation with no tools required

## TECHNICAL CHARACTERISTICS

**FILT'RAY 1-month tap filter - Ref. AT030728**

|            |                                                                                                                                                                                    |
|------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Technology | Point-of-use filtration                                                                                                                                                            |
| Height     | 72mm                                                                                                                                                                               |
| Length     | 75mm                                                                                                                                                                               |
| Flow rate  | 6 lpm at 3 bar, 4 lpm at 1 bar                                                                                                                                                     |
| Lifespan   | 1 month                                                                                                                                                                            |
| Norms      | <b>ACS</b>   |

