

# FILT'RAY Compact 1-month shower head filter

Ref. AT30161

Non-sterile

## DESCRIPTION

**FILT'RAY Compact 1-month shower head filter**

**Ref. AT30161**

AQUATOOLS anti-Legionella and all-germ shower filter.  
Device made from ABS.  
Membrane type: tubular microfiltration.  
Membrane made from polyethylene.  
Pore cut-off threshold: 0.1µm nominal and 0.2µm absolute sterilising grade.  
Qualified water volume per duration of use: 3,000 litres.  
Filtration flow rate: 6 lpm at 1 bar, 11.2 lpm at 3 bar and 13 lpm at 5 bar (filter only, excludes any flow rate restriction in the shower valve).  
Maximum upstream pressure at the point-of-use: 5 bar.  
Maximum lifespan: 1 month after installation.  
Resistance to thermal shocks: 70°C for 30 minutes.  
Resistance to chemical shocks: active chlorine levels up to 100 ppm for 1 hour (cumulative total).  
Non-sterile filter.  
Medical device CE Class I.  
Weight: empty filter 105g / full filter 186g.  
Reference for a single unit. Must be ordered in multiples of 10.  
Packaged in boxes of 10 units.

## TECHNICAL CHARACTERISTICS

**FILT'RAY Compact 1-month shower head filter - Ref. AT30161**

|            |                                                       |
|------------|-------------------------------------------------------|
| Connector  | M1/2"                                                 |
| Technology | Point-of-use filtration                               |
| Height     | 62mm                                                  |
| Length     | 194mm                                                 |
| Flow rate  | 6 lpm at 1 bar, 11.2 lpm at 3 bar and 13 lpm at 5 bar |
| Diameter   | 56mm                                                  |
| Lifespan   | 1 month                                               |

Norms **ACS** **CE** **WRAS** **KTW BWGL**



## ADVANTAGES



Anti-Legionella and all-germ



Space saving



Lifespan: 1 month after installation



Quick to install without tools

